


Dura-line XD 17 and XD⁺ 17

Robust, cost-efficient, wireless flat panel detector

Versatile, robust and cost-effective, Dura-line XD 17 and XD⁺ 17 wireless DR detectors provide exquisite image quality, with a 17 x 17" (43 x 43 cm) field of view, choice of 1 or 2 batteries for extended battery autonomy, and a Cesium Iodide scintillator.



DURA-LINE FLAT PANEL DETECTOR

BATTERY PACK

- Model: FXRB-04A
- Type: Lithium Ion Polymer
- Nominal Voltage: DC +11.55V
- Nominal Capacity: 3,400mAh
- Number of Cells: 3S1P (3 Series 1 Parallel)
- Battery Life: Approx. 800 times
(One time standard: Fully charged and then fully discharged)
- Dimensions (H × W × D):
189.0mm × 89.0mm × 6.65mm (max.)
- Weight: 185g (max.)

BATTERY CHARGER, CAN CHARGE UP TO 2 DETECTOR BATTERIES


- Model: FXRC-04A
- Nominal Voltage: DC +24V
- Dimensions (H × W × D): 230.0mm × 304.0mm × 15.5mm
- Weight: 450gr

CRADLE CHARGER, CAN CHARGE UP TO 2 DETECTORS

- Model: FXRR-01A
- Dimensions (H × W × D): 159.0mm × 410.0mm × 168.0mm
- Weight: 4.04kg



DQE SPECIFICATIONS:



DURA-LINE FLAT PANEL DETECTOR

ORDER CODES (ABC CODES)

■ MAXY0	XD 17 DETECTOR (FXRD-4343VAW)
■ MAX2U	XD+ 17 DETECTOR (FXRD-4343VAW PLUS)
■ MA3YA	BATTERY CHARGER DURALINE (FXRC-04A)
■ MAX87	BATTERY PACK DURALINE (FXRB-04A)
■ MAX4Y	CRADLE CHARGER DURALINE (FXRR-01A)
■ MAX63	MINI-SCU BOX DURALINE (FXRS-04A)
■ MAX51	SCU-LITE DURALINE (FXRP-02A)
■ MAX75	TETHERED CABLE XD
■ MA3X7	XD/XD+ 17 CLICK ON HANDLE
■ MA3V3	XD/XD+ 17 GRID R 10:1, 132 L/INCH
■ MA3U1	XD/XD+ 17 GRID, R 6:1, 132 L/INCH

For more information on Agfa, please visit our website on www.agfa.com ■

Agfa and the Agfa rhombus are trademarks of Agfa-Gevaert NV, Belgium, or its affiliates. MUSICA is a trademark of Agfa-Gevaert NV, Belgium, or its affiliates. All rights reserved. All information contained herein is intended for guidance purposes only, and characteristics of the products and services described in this publication can be changed at any time without notice. Products and services may not be available for your local area. Please contact your local sales representative for availability information. Agfa-Gevaert NV diligently strives to provide as accurate information as possible, but shall not be responsible for any typographical error.

© 2021 Agfa NV
All rights reserved
Published by Agfa NV
Septestraat 27 - 2640 Mortsel
Belgium

GB 00202105

AGFA 