



# Acuity

## FDR+



**WORKFLOW INNOVATION**

Radmedix 2510 Lance Drive, Kettering, OH 45409 USA  
Tel: 844.723.6334 | Fax: 513.278.0150 | [info@radmedix.com](mailto:info@radmedix.com)  
© 2021 Radmedix Inc. All rights reserved.



# FDR+ ADVANCED SYSTEM

Acuity FDR+ x-ray system with a standard configuration of floor mounted tube stand, vertical bucky stand and bucky table. This system has advanced safety workflow and positioning features. It also has an intuitive, 10" user-friendly touchscreen that provides seamless operation of your radiographic system. It displays graphical tube angulation, SID values, and tube to bucky auto sync. With such a wide range of features, the Acuity FDR+ is considered our highest performing x-ray system offering unparalleled price-performance ratio.

## 10" Intuitive Touchscreen

Preset source to image distance with Auto Sync feature.

Real time display of SID and tube angulation.

Suitable for any clinical application.



## Premium Wall Stand

Motorized and counterbalance adjustment controls with electromagnetic locks.

Convenient handgrips to optimize patient positioning.

Overhead handgrips to promote proper upright exams.



AEC Optional

# FDR+ KEY HIGHLIGHTS



Auto Sync Tracking with Wall & Table Bucky's



Rotational Column



6 Way Elevating Float Top Table



Wireless Pedal



AEC Optional

Emergency Stop

# EFFORTLESS POSITIONING

Motorized auto-tracking enables a wide array of patient positioning.



# INTUITIVE 10" TOUCHSCREEN



Elevated Table Specifications  
Lowering Height: 21 in  
Patient Load: 550 lbs  
Tabletop dimensions: 32 x 86 1/2 in  
Integrated Emergency Stop

# FDR STANDARD SYSTEM

The Acuity FDR is also offered in our standard configuration. Offering the same versatility of the Acuity FDR+, but with our standard non-touchscreen display and advanced features. The Acuity FDR Standard offers the same exceptional digital imaging quality but with standard features at a more budget friendly price point.

## 4" Console Display

Preset source to image distance with Auto Sync feature.

Real time display of SID and tube angulation.

Suitable for any clinical application.



AEC Optional

## Premium Wall Stand

Manual counterbalance adjustment controls with electromagnetic locks.

Convenient handgrips to optimize patient positioning

Overhead handgrips to promote proper upright exams.



# FDR STANDARD FEATURES



Rotational Column



4 Way Float Non Elevating Table



Built In Foot Pedal



AEC Optional